

Freedom® Total Knee System



Carefully read all instructions and be familiar with the surgical techniques prior to use.

Description

The Maxx Orthopedics' Freedom® Total Knee System is a system of components intended to replace the femoral, tibial and patella articular surfaces of the knee joint. Components are available in many styles and sizes and are manufactured from various types of metals and non-metallic materials. The different product categories include:

1. The Freedom Total Knee System
 - a. Metal Backed Tibial Components
 - b. All-Poly Tibial Components (CR and PS)
 - c. Femoral (CR and PS)
2. Freedom Stemmed Tibial Components
3. Freedom PCK Components

The component style, size, compatibility, and specific component material is provided on the outside carton label. All implantable components are designed for single use only. Note: PCK components are not licensed for sale in Canada.

Product Selection Information

- Appropriate matching of the components will occur when the articular component is matched to the femoral component (by letter size designation and constraint style) and to the tibial base (by numerical size designation). For example, a size **B1-2 PS (Posterior Stabilized)** articular component is appropriately matched to a size **B PS** femoral component and either a size **1** or size **2** tibial base plate. The UC (Ultra Congruent) constraint style tibial articular components are matched with corresponding letter size CR femoral components. Mismatching may result in poor surface contact and produce pain, decrease wear resistance, produce instability of the implant, or otherwise reduce implant life.
- Use only instruments and trials specifically designed for use with these devices to help ensure accurate surgical implementation, soft tissue balancing, and evaluation of knee function.
- Selections between the various sizes is a matter of physician discretion.
- Do not use products past their expiration date.

Indications

- Severe knee joint pain, loss of mobility, and disability due to: rheumatoid arthritis, osteoarthritis, traumatic arthritis, polyarthritis.
- Correction of functional deformities.
- Post-traumatic loss of knee joint contour, particularly when there is patellofemoral erosion, dysfunction, or prior patellectomy.
- Moderate valgus, varus, or flexion trauma.
- Knee fractures untreatable by other methods.
- Revision surgery where sufficient bone stock and soft tissue integrity are present is indicated only for PCK components and Primary PCK components.

The Freedom Total Knee System, Freedom Stemmed Tibial Components and Freedom PCK Components are indicated for cemented fixation.

Contraindications

- History of infection in the affected joint that may affect the function of the implanted prosthetic.
- Less than optimal bone stock on femoral or tibial surfaces resulting from a history of disease, infection, or prior surgical procedures which cannot provide adequate support for the implantation.
- Compromised skeletal bone quality.
- Neuropathic disease that adversely affects the prosthetic joint.
- Osteoporosis or deficiency of musculature that compromises the affected limb.
- Pain-free and stable arthrodesis in an adequate functional position.
- Unstable knee joint secondary to negative collateral ligament integrity. The conditions like obesity or overweight*, active sports participation, high levels of patient activity**, tend to impose severe loading on the affected extremity, thereby placing the patient at higher risk for failure of the knee implant.
- The safety and effectiveness of the use of the Freedom total knee system in patients depends on various factors which include, but are not limited to, surgical technique, patient build, preoperative flexion and age.
- Known allergic reactions to metallic corrosions, wear debris, or cement particles.

Note:

* WHO (World Health Organization) defines "overweight" as a BMI (Body mass index) greater than or equal to 25, and "obesity" as a BMI greater than or equal to 30.

**According to The Knee Society, appropriate activities for patients following joint replacement surgery include cycling, calisthenics, swimming, low-resistance rowing, walking, hiking, low-resistance weightlifting and use of stationary skiing machines. Other suitable activities include bowling, croquet, golf, doubles tennis, table tennis and ballroom and square dancing. According to The Knee Society, activities such as downhill skiing, scuba diving, in-line skating, ice skating, softball, volleyball, speed walking, horseback riding, hunting and low-impact aerobics are suitable, but include more risk. Baseball, basketball, football, hockey, soccer, gymnastics, jogging, rock climbing, hang gliding, parachuting and high-impact aerobics are activities that should be avoided.

^a Ibrahim Akkawi et al. Sports after Total Knee Prosthesis. Springer-Verlag Berlin Heidelberg 2014.

Adverse Effects

- Long term swelling or infection.
- No improvement in range of motion.
- Neuropathic disorders.
- Dislocations, bone fractures, and/or joint instability.
- Wear and/or deformation of articulating surfaces.
- Venous thrombosis.
- Prolonged and excessive joint pain and/or inflammation.

- Aseptic loosening of implant.
- Possible detachment of the coating(s) on components with porous coating, potentially leading to increased debris particles

Warnings and Precautions:

- This device is intended for cemented use only. Whereby stem extensions are either for cemented or press-fit use. Stems with length of 40 mm are recommended to be used cemented.
- This device is for single patient use only. Never reuse the implant, even though it may appear undamaged as the implant may have developed microscopic imperfections which could lead to failure. Reuse can potentially compromise device performance and patient safety. If prosthesis is reused, there are chances of infection, loosening or revision surgery.
- Malalignment of the components or inaccurate implantation can lead to excessive wear and/or failure of the implant or procedure.
- Polished bearing areas must not come in contact with hard or abrasive surfaces. Bearing areas must be free of debris and clean prior to assembly.
- Contouring or bending of the implant may reduce its fatigue strength and cause premature failure under load.
- All-polyethylene tibial base plates should be limited to use in low demand (i.e. modest weight, activity level) patients with good bone quality.
- Return all packages with flaws in the sterile barrier to the supplier. Do not re-sterilize.
- The optional Freedom Stemmed Tibial Augments are intended for screw attachment to the Stemmed Tibial Base Plate.
- The optional Freedom PCK Femoral Augments are intended for screw attachment to the PCK Femoral Component.
- The Freedom branded devices have not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Freedom branded devices in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
- In case of revision surgery, special care should be taken on removal of primary device and while implanting revision device. When removing the components to be revised, care must be taken to preserve as much of the remaining bone stock as possible and to avoid the risk of fracture of the residual bone.

Storage Condition:

To be stored in temperature range of 15 °C to 25 °C.

Materials

All of the materials used in the Maxx Orthopedics' Freedom® Total Knee System meet ASTM standards for implants.

Sterilization

- All Implants are supplied sterile in protective packaging to a Sterility Assurance Level (SAL) of 10⁻⁶ using either gamma radiation or Ethylene Oxide.
- The method of sterilization is provided on the outside carton label. Components sterilized through radiation have been exposed to a minimum dose of 25 kGy of gamma radiation.
- The packaging of all sterile products should be inspected for flaws in the sterile barrier before opening. In the event of such a flaw, the product must be returned non-sterile. Trial components should be used to avoid having to open any aspect of the sterile package prior to the component use.
- Care should be taken to prevent contamination of the component. In the event of contamination, this product should be discarded.
- If the package is opened, but the product is not used, the component should not be re-sterilized and should be discarded or returned to the supplier.

Component Type	Material	Applicable standard	Sterilization Method
Femoral Components (CR, PS and PCK)	Cobalt Chromium Molybdenum (CoCrMo) alloy	ISO 5832-4/ ASTM F75	Ethylene Oxide
Tibial Articular Surfaces (CR, PS, and PCK)	Ultra-High Molecular Weight Polyethylene (UHMWPE)	ISO 5834-4/ ASTM F648	Ethylene Oxide
Tibial Base Plates (Standard and Stemmed)	Cobalt Chromium Molybdenum (CoCrMo) alloy	ISO 5832-4 / ASTM F75	Ethylene Oxide
All Poly Patellas	Ultra-High Molecular Weight Polyethylene (UHMWPE)	ISO 5834-4/ ASTM F648	Ethylene Oxide
Offset Junctions	Ti-6Al-4V-ELI	ISO 5832-3/ ASTM F136	Gamma Irradiation
Stem Extensions	Ti-6Al-4V-ELI	ISO 5832-3/ ASTM F136	Gamma Irradiation
Augments (Tibial and Femoral)	Ti-6Al-4V-ELI	ISO 5832-3/ ASTM F136	Gamma Irradiation

Patient Guidance

The probability of resulting complications and/or failure of total knee prostheses is increased in cases where the patient's physical and medical presentation (e.g. weight, or other diseases) is a detriment, the patient's functional outcome goals are above what is reasonably attainable, the patient's expectations for knee joint function are unrealistic, and/or the patient's engages in a less than optimal level of rehabilitation post-surgery. Excessive physical activity should be limited. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail, or need to be replaced. The implant may not, and is not guaranteed to last the rest of the patient's life. Because prosthetic joints are not as strong, reliable, or durable as natural, healthy joints, all prosthetic knees may need to be replaced at some point. The prosthesis will not restore functions to the level expected with normal healthy bone, and the patient should temper their expectations to a realistic level. As with all prosthetic implants, the durability of the components is affected by numerous biologic, biomechanical, and other external factors which can limit their service life. Adherence to the indications, contraindications, precautions, and warnings for this product is essential for maximizing service life.

CAUTION: FEDERAL LAW (U.S.A) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN.

Freedom® Total Knee System

R_x



Maxx Orthopedics, Inc.
2460 General Armistead Avenue
Suite 100, Norristown, PA
19403, USA
Tel: +1-484-342-0092
Fax: +1-484-342-0096
Email: info@maxxortho.com

CE
2460

EC REP

AJW Consulting GmbH, [Breite Straße 3](https://www.breite-strasse.com), 40213 Dusseldorf, Germany
Phone: +49) 211 54059 6030
E-mail: info@ajwtech.com

Carefully read all instructions and be familiar with the surgical techniques prior to use.

SYMBOL KEY LEGEND	
	MANUFACTURER
REF	CATALOG NUMBER
LOT	BATCH CODE
	DATE OF MANUFACTURE
	USE BY DATE
	CAUTION
	DO NOT RE-USE
STERILE EO	STERILIZED WITH ETYLENE OXIDE
STERILE R	STERILIZED WITH GAMMA IRRADIATION
EC REP	EUROPEAN AUTHORISED REPRESENTATIVE
R _x	PRESCRIPTION ONLY
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY
	DO NOT USE IF PACKAGE IS DAMAGE
	TEMPERATURE LIMITATION