


maxx

orthopedics



Freedom® Medial Congruent Insert

System Overview


F R E E D O M[®]



K N E E

Freedom[®] Medial Congruent Insert


The Freedom Medial Congruent Insert optimizes knee stability and replicates native knee function in total knee arthroplasty. The asymmetric condylar design enables lateral femoral rollback, while the increased anterior lip prevents subluxation in the medial compartment.



- Medial congruency for better stability throughout full ROM
- Enhanced stability in extension and greater translation in lateral compartment.
- Designed to replicate normal knee kinematics.
- Increased anterior lip prevents femoral subluxation with or without PCL.
- Compatible with cruciate retaining and cruciate sacrificing techniques.
- UHMWPE GUR 1020 for excellent wear resistance and durability
- Available in 9, 10, 11, 12, 14, & 17mm.


Device Sizing

		Femoral Component Size							
		A	B	C	D	E	F	G	H
Tibial Baseplate Size	1	1-2	1-2	1-2	1-2				
	2								
	3			3-4	3-4	3-4	3-4		
	4								
	5					5-6	5-6	5-6	
	6								
	7						7-8	7-8	7-8
	8								




Tibial Insert Size	Anterior Medial Lip Height (mm)	Posterior Medial Lip Height (mm)
A1-2	9.0	3.0
B1-2	9.5	3.0
C1-2	10.0	3.0
C3-4	10.0	3.0
D1-2	10.5	3.1
D3-4	10.5	3.1
E3-4	11.0	3.2
E5-6	11.0	3.2
F3-4	11.5	3.2
F5-6	11.5	3.2
F7-8	11.5	3.2
G5-6	12.0	3.3
G7-8	12.0	3.3
H7-8	12.5	3.4

Set Configuration




INS-PK-MCF




Insert & Femoral Trials

INS-PK-MCF (Top Tray)




Part #	Description
CBLA00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE A, MC, 9mm Compatible, LEFT
CBRA00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE A, MC, 9mm Compatible, RIGHT
CBLB00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE B, MC, 9mm Compatible, LEFT
CBRB00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE B, MC, 9mm Compatible, RIGHT
CBLC00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE C, MC, 9mm Compatible, LEFT
CBRC00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE C, MC, 9mm Compatible, RIGHT
CBLD00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE D, MC, 9mm Compatible, LEFT
CBRD00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE D, MC, 9mm Compatible, RIGHT
CBLE00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE E, MC, 9mm Compatible, LEFT
CBRE00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE E, MC, 9mm Compatible, RIGHT
CBLF00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE F, MC, 9mm Compatible, LEFT
CBRF00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE F, MC, 9mm Compatible, RIGHT
CBLG00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE G, MC, 9mm Compatible, LEFT
CBRG00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE G, MC, 9mm Compatible, RIGHT
CBLH00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE H, MC, 9mm Compatible, LEFT
CBRH00MC-F	TIBIAL TRIAL INSERT (TOP), SIZE H, MC, 9mm Compatible, RIGHT
CBX109MC-F	TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 9mm
CBX110MC-F	TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 10mm
CBX111MC-F	TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 11mm
CBX112MC-F	TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 12mm
CBX114MC-F	TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 14mm
CBX117MC-F	TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 17mm
CBX309MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 9mm
CBX310MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 10mm

INS-PK-MCF (Top Tray) - Continued



Part #	Description
CBX311MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 11mm
CBX312MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 12mm
CBX314MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 14mm
CBX317MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 17mm
CBX509MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 9mm
CBX510MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 10mm
CBX511MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 11mm
CBX512MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 12mm
CBX514MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 14mm
CBX517MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 17mm
CBX709MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 9mm
CBX710MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 10mm
CBX711MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 11mm
CBX712MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 12mm
CBX714MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 14mm
CBX717MC-F	TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 17mm

INS-PK-MCF (Bottom Tray)



Part #	Description
CCLB01CR-K	PEGLESS FEMORAL TRIAL SIZE B, CR, L
CCLC01CR-K	PEGLESS FEMORAL TRIAL SIZE C, CR, L
CCLD01CR-K	PEGLESS FEMORAL TRIAL SIZE D, CR, L
CCLE01CR-K	PEGLESS FEMORAL TRIAL SIZE E, CR, L
CCLF01CR-K	PEGLESS FEMORAL TRIAL SIZE F, CR, L
CCLG01CR-K	PEGLESS FEMORAL TRIAL SIZE G, CR, L
CCLH01CR-K	PEGLESS FEMORAL TRIAL SIZE H, CR, L
CCRB01CR-K	PEGLESS FEMORAL TRIAL SIZE B, CR, R
CCRC01CR-K	PEGLESS FEMORAL TRIAL SIZE C, CR, R
CCRD01CR-K	PEGLESS FEMORAL TRIAL SIZE D, CR, R
CCRE01CR-K	PEGLESS FEMORAL TRIAL SIZE E, CR, R
CCRF01CR-K	PEGLESS FEMORAL TRIAL SIZE F, CR, R
CCRG01CR-K	PEGLESS FEMORAL TRIAL SIZE G, CR, R
CCRH01CR-K	PEGLESS FEMORAL TRIAL SIZE H, CR, R

Pursue Life

For more information about Freedom® Knee System please contact your local representative.

Find us on social media



@maxxortho

Download the Maxx Ortho App for more product information:



Download on the
App Store



GET IT ON
Google Play

Search for Maxx Ortho

Freedom® Total Knee System



Maxx Orthopedics, Inc.
2460 General Armistead Ave, Ste 100
Norristown, Pa 19403 USA



Carefully read all instructions and be familiar with the surgical techniques prior to use.

Please see the package insert for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with Freedom Total Knee® System

CAUTION: THIS DEVICE IS RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

Freedom Knee® is manufactured by Maxx Orthopedics, Inc. FREEDOM, FREEDOM KNEE, PCK, RENEW, Pursue Life and Think Outside the Box are Registered Trademarks of Maxx Orthopedics, Inc.

©2025 Maxx Orthopedics. All rights reserved.

MXO-MP00058 R01 03/25

maxx
orthopedics

www.maxxortho.com