


FREEDOM® MEDIAL CONGRUENT INSERT

The Freedom Medial Congruent Insert is designed to optimize knee stability and replicate native knee function in total knee arthroplasty. The asymmetric condylar design enables lateral femoral rollback, while the increased anterior lip prevents subluxation in the medial compartment.




- Medial congruency for better stability throughout full ROM
- Enhanced stability in extension and greater translation in lateral compartment.
- Designed to replicate normal knee kinematics.
- Increased anterior lip prevents femoral subluxation with or without PCL.
- Compatible with cruciate retaining and cruciate sacrificing techniques.



“Mimics native knee kinematics”
- Dr. Siraj Sayeed

Device Sizing

| | | Femoral Component Size | | | | | | | |
|-----------------------|---|------------------------|-----|-----|-----|-----|-----|-----|-----|
| | | A | B | C | D | E | F | G | H |
| Tibial Baseplate Size | 1 | 1-2 | 1-2 | 1-2 | 1-2 | | | | |
| | 2 | | | | | | | | |
| | 3 | | | 3-4 | 3-4 | 3-4 | 3-4 | | |
| | 4 | | | | | | | | |
| | 5 | | | | | 5-6 | 5-6 | 5-6 | |
| | 6 | | | | | | | | |
| | 7 | | | | | | 7-8 | 7-8 | 7-8 |
| | 8 | | | | | | | | |




| Tibial Insert Size | Anterior Medial Lip Height (mm) | Posterior Medial Lip Height (mm) |
|--------------------|---------------------------------|----------------------------------|
| A1-2 | 9.0 | 3.0 |
| B1-2 | 9.5 | 3.0 |
| C1-2 | 10.0 | 3.0 |
| C3-4 | 10.0 | 3.0 |
| D1-2 | 10.5 | 3.1 |
| D3-4 | 10.5 | 3.1 |
| E3-4 | 11.0 | 3.2 |
| E5-6 | 11.0 | 3.2 |
| F3-4 | 11.5 | 3.2 |
| F5-6 | 11.5 | 3.2 |
| F7-8 | 11.5 | 3.2 |
| G5-6 | 12.0 | 3.3 |
| G7-8 | 12.0 | 3.3 |
| H7-8 | 12.5 | 3.4 |

Available in 9, 10, 11, 12, 14, & 17mm insert thickness.

Set Configuration



INS-PK-MCF (Trial Inserts)



INS-PK-MCF (Femoral Trials)

INS-PK-MCF (Top Tray)




| Part # | Description |
|------------|--|
| CBLA00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE A, MC, 9mm Compatible, LEFT |
| CBRA00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE A, MC, 9mm Compatible, RIGHT |
| CBLB00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE B, MC, 9mm Compatible, LEFT |
| CBRB00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE B, MC, 9mm Compatible, RIGHT |
| CBLC00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE C, MC, 9mm Compatible, LEFT |
| CBRC00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE C, MC, 9mm Compatible, RIGHT |
| CBLD00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE D, MC, 9mm Compatible, LEFT |
| CBRD00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE D, MC, 9mm Compatible, RIGHT |
| CBLE00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE E, MC, 9mm Compatible, LEFT |
| CBRE00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE E, MC, 9mm Compatible, RIGHT |
| CBLF00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE F, MC, 9mm Compatible, LEFT |
| CBRF00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE F, MC, 9mm Compatible, RIGHT |
| CBLG00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE G, MC, 9mm Compatible, LEFT |
| CBRG00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE G, MC, 9mm Compatible, RIGHT |
| CBLH00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE H, MC, 9mm Compatible, LEFT |
| CBRH00MC-F | TIBIAL TRIAL INSERT (TOP), SIZE H, MC, 9mm Compatible, RIGHT |
| CBX109MC-F | TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 9mm |
| CBX110MC-F | TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 10mm |
| CBX111MC-F | TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 11mm |
| CBX112MC-F | TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 12mm |
| CBX114MC-F | TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 14mm |
| CBX117MC-F | TIBIAL TRIAL INSERT BOTTOM, SIZE 1-2, 9mm COMPATIBLE, 17mm |
| CBX309MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 9mm |
| CBX310MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 10mm |

INS-PK-MCF (Top Tray) - Continued



| Part # | Description |
|------------|--|
| CBX311MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 11mm |
| CBX312MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 12mm |
| CBX314MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 14mm |
| CBX317MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 3-4, 9mm COMPATIBLE, 17mm |
| CBX509MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 9mm |
| CBX510MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 10mm |
| CBX511MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 11mm |
| CBX512MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 12mm |
| CBX514MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 14mm |
| CBX517MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 5-6, 9mm COMPATIBLE, 17mm |
| CBX709MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 9mm |
| CBX710MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 10mm |
| CBX711MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 11mm |
| CBX712MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 12mm |
| CBX714MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 14mm |
| CBX717MC-F | TIBIAL TRIAL INSERT BOTTOM, Size 7-8, 9mm COMPATIBLE, 17mm |

INS-PK-MCF (Bottom Tray)



| Part # | Description |
|------------|-------------------------------------|
| CCLB01CR-K | PEGLESS FEMORAL TRIAL SIZE B, CR, L |
| CCLC01CR-K | PEGLESS FEMORAL TRIAL SIZE C, CR, L |
| CCLD01CR-K | PEGLESS FEMORAL TRIAL SIZE D, CR, L |
| CCLE01CR-K | PEGLESS FEMORAL TRIAL SIZE E, CR, L |
| CCLF01CR-K | PEGLESS FEMORAL TRIAL SIZE F, CR, L |
| CCLG01CR-K | PEGLESS FEMORAL TRIAL SIZE G, CR, L |
| CCLH01CR-K | PEGLESS FEMORAL TRIAL SIZE H, CR, L |
| CCRB01CR-K | PEGLESS FEMORAL TRIAL SIZE B, CR, R |
| CCRC01CR-K | PEGLESS FEMORAL TRIAL SIZE C, CR, R |
| CCRD01CR-K | PEGLESS FEMORAL TRIAL SIZE D, CR, R |
| CCRE01CR-K | PEGLESS FEMORAL TRIAL SIZE E, CR, R |
| CCRF01CR-K | PEGLESS FEMORAL TRIAL SIZE F, CR, R |
| CCRG01CR-K | PEGLESS FEMORAL TRIAL SIZE G, CR, R |
| CCRH01CR-K | PEGLESS FEMORAL TRIAL SIZE H, CR, R |

